
Product Manual


Cdc42 Activation Assay Kit

Catalog Number

STA-402

20 assays

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures



CELL BIOLABS, INC.
Creating Solutions for Life Science Research

Introduction

Small GTP-binding proteins (or GTPases) are a family of proteins that serve as molecular regulators in signaling transduction pathways. Cdc42, a 21 kDa protein, belongs to the family of Rho GTPases regulates a variety of biological response pathways that include cell motility, cell division, gene transcription, and cell transformation. Like other small GTPases, Cdc42 regulates molecular events by cycling between an inactive GDP-bound form and an active GTP-bound form. In its active (GTP-bound) state, Cdc42 binds specifically to the p21-binding domain (PBD) of p21-activated protein kinase (PAK) to control downstream signaling cascades.


Cell Biolabs' Cdc42 Activation Assay Kit utilizes PAK PBD Agarose beads to selectively isolate and pull-down the active form of Cdc42 from purified samples or endogenous lysates. Subsequently, the precipitated GTP-Cdc42 is detected by western blot analysis using an anti-Cdc42 specific monoclonal antibody (see Figure 3 and Assay Principle).

Cell Biolabs' Cdc42 Activation Assay Kit provides a simple and fast tool to monitor the activation of Cdc42. The kit includes easily identifiable PAK1 PBD Agarose beads (see Figure 1), pink in color, and a Cdc42 Immunoblot Positive Control for quick Cdc42 identification. Each kit provides sufficient quantities to perform 20 assays.



Figure 1: PAK RBD Agarose beads, in color, are easy to visualize, minimizing potential loss during washes and aspirations.

Assay Principle



Related Products

1. STA-400: Pan-Ras Activation Assay Kit
2. STA-401-1: Rac1 Activation Assay
3. STA-401-2: Rac2 Activation Assay
4. STA-403-A: RhoA Activation Assay Kit
5. STA-405: RhoA/Rac1/Cdc42 Activation Assay Combo Kit

Kit Components (shipped on blue ice)

1. PAK1 PBD Agarose (Part No. STA-411): One vial – 800 μ L of 50% slurry, 400 μ g PBD in PBS containing 50% glycerol.
Note: Agarose bead appears pink in color for easy identification, washing, and aspiration.
2. 100X GTP γ S (Part No. 240103): One vial – 50 μ L of 10 mM GTP γ S dissolved in sterile water.
3. 100X GDP (Part No. 240104): One vial – 50 μ L of 100 mM GDP dissolved in sterile water.
4. 5X Assay/Lysis Buffer (Part No. 240102): One bottle – 30 mL of 125 mM HEPES, pH 7.5, 750 mM NaCl, 5% Igepal CA-630, 50 mM MgCl₂, 5 mM EDTA, 10% Glycerol.
5. Anti-Cdc42, Mouse Monoclonal (Part No. 240201): One vial – 40 μ L in PBS, pH 7.4, 0.05% NaN₃, 0.1% BSA.
Note: This monoclonal antibody specifically reacts with human, mouse, and rat Cdc42.
6. Cdc42 Immunoblot Positive Control (Part No. 240210): One vial – 100 μ L of partially purified, recombinant Cdc42 from *E. coli* (provided ready-to-use in 1X reducing SDS-PAGE Sample Buffer, pre-boiled).

Materials Not Supplied

1. Stimulated and non-stimulated cell lysates
2. Cdc42 activators
3. Protease inhibitors
4. 0.5 M EDTA in water
5. 1 M MgCl₂
6. 30°C incubator or water bath
7. 4°C tube rocker or shaker
8. 2X reducing SDS-PAGE sample buffer
9. Electrophoresis and immunoblotting systems
10. Immunoblotting wash buffer such as TBST (10 mM Tris-HCl, pH 7.4, 0.15 M NaCl, 0.05% Tween-20)
11. Immunoblotting blocking buffer (TBST containing 5% Non-fat Dry Milk)
12. PVDF or nitrocellulose membrane
13. Secondary Antibody
14. ECL Detection Reagents

Storage

Store all kit components at -20°C. The 5X Assay/Lysis Buffer may be stored at either -20°C or 4°C. Avoid multiple freeze/thaw cycles.

Preparation of Reagents

- 1X Assay/Lysis Buffer: Mix the 5X Stock briefly and dilute to 1X in deionized water. Just prior to usage, add protease inhibitors such as 1 mM PMSF, 10 µg/mL leupeptin, and 10 µg/mL aprotinin.

Preparation of Samples

Note: It is advisable to use fresh cell lysates because GTP-Cdc42 is quickly hydrolyzed to GDP-Cdc42; frozen lysates stored at -70°C may be used. Performing steps at 4°C or on ice may reduce hydrolysis. Avoid multiple freeze/thaw cycles of lysates.

I. Adherent Cells

1. Culture cells to approximately 80-90% confluence. Stimulate cells with Cdc42 activator(s) as desired.
2. Aspirate the culture media and wash twice with ice-cold PBS.
3. Completely remove the final PBS wash and add ice-cold 1X Assay/Lysis Buffer to the cells (0.5 - 1 mL per 100 mm tissue culture plate).
4. Place the culture plates on ice for 10-20 minutes.
5. Detach the cells from the plates by scraping with a cell scraper.
6. Transfer the lysates to appropriate size tubes and place on ice.
7. If nuclear lysis occurs, the cell lysates may become very viscous and difficult to pipette. If this occurs, lysates can be passed through a 27½-gauge syringe needle 3-4 times to shear the genomic DNA.
8. Clear the lysates by centrifugation for 10 minutes (14,000 x g at 4°C).
9. Collect the supernatant and store samples on ice for immediate use, or snap freeze and store at -70°C for future use.
10. Proceed to GTPγS/GDP Loading for positive and negative controls, or Pull-Down Assay.

II. Suspension Cells

1. Culture cells and stimulate with Cdc42 activator(s) as desired.
2. Perform a cell count, and then pellet the cells by centrifugation.
3. Aspirate the culture media and wash twice with ice-cold PBS.
4. Completely remove the final PBS wash and add ice-cold 1X Assay/Lysis Buffer to the cell pellet (0.5 – 1 mL per 1×10^7 cells).
5. Lyse the cells by repeated pipetting.
6. Transfer the lysates to appropriate size tubes and place on ice.
7. If nuclear lysis occurs, the cell lysates may become very viscous and difficult to pipette. If this occurs, lysates can be passed through a 27½-gauge syringe needle 3-4 times to shear the genomic DNA.

8. Clear the lysates by centrifugation for 10 minutes (14,000 x g at 4 °C).
9. Collect the supernatant and store samples on ice for immediate use, or snap freeze and store at -70 °C for future use.
10. Proceed to GTP γ S/GDP Loading for positive and negative controls, or Pull-Down Assay.

Assay Protocol

Important Note: Before running any Small GTPase pulldown assay, it is always a good practice to run a Western Blot directly on the cell lysate using the antibody provided in this kit. For example: load 5 μ g, 10 μ g and 20 μ g of lysate onto an SDS-PAGE gel, transfer and blot. When proceeding with the pulldown assay, use 100-times the amount of lysate that gave you a clear band of your desired small GTPase in the direct Western blot. For example: if the 5- μ g band was faint but the 10- μ g band was clear and strong, use 100 x 10 μ g = 1 mg of lysate in the assay. Using sufficient lysate in the pulldown assay is critical to success.

I. GTP γ S/GDP Loading (Positive and Negative Controls)

Note: Samples that will not be GTP γ S/GDP loaded may be kept on ice during the loading of controls.

1. Aliquot 0.5 – 1 mL of each cell lysate to two microcentrifuge tubes.
Note: Typical protein content/sample is > 0.5 mg.
2. Adjust the volume of each sample to 1 mL with 1X Assay Lysis Buffer.
3. Add 20 μ L of 0.5 M EDTA to each sample.
4. Add 10 μ L of 100X GTP γ S to one tube (positive control) and 10 μ L of 100X GDP to the other tube (negative control). Mix and label each tube appropriately.
5. Incubate the tubes for 30 minutes at 30°C with agitation.
6. Stop the loading by adding 65 μ L of 1 M MgCl₂ to each tube. Mix and place tubes on ice.
7. Continue with Pull-Down assay.

II. Cdc42 Pull-Down Assay

1. Aliquot 0.5 – 1 mL of cell lysate (treated with Cdc42 activators or untreated) to a microcentrifuge tube.
2. Adjust the volume of each sample to 1 mL with 1X Assay Lysis Buffer.
3. Thoroughly resuspend the PAK PBD Agarose bead slurry by vortexing or titrating.
4. Quickly add 40 μ L of resuspended bead slurry to each tube (including GTP γ S/GDP controls).
5. Incubate the tubes at 4°C for 1 hour with gentle agitation.
6. Pellet the beads by centrifugation for 10 seconds at 14,000 x g.
7. Aspirate and discard the supernatant, making sure not to disturb/remove the bead pellet.
8. Wash the bead 3 times with 0.5 mL of 1X Assay Buffer, centrifuging and aspirating each time.
9. After the last wash, pellet the beads and carefully remove all the supernatant.
10. Resuspend the bead pellet in 40 μ L of 2X reducing SDS-PAGE sample buffer.

11. Boil each sample for 5 minutes.
12. Centrifuge each sample for 10 seconds at 14,000 x g.

III. Electrophoresis and Transfer

1. Load 20 μ L/well of pull-down supernatant to a polyacrylamide gel. Also, it's recommended to include a pre-stained MW standard (as an indicator of a successful transfer in step 3).

Note: If desired, 10 μ L/well of Cdc42 Immunoblot Positive Control (provided ready-to-use, pre-boiled) can be added as an immunoblot positive control.

2. Perform SDS-PAGE as per the manufacturer's instructions.
3. Transfer the gel proteins to a PVDF or nitrocellulose membrane as per the manufacturer's instructions.

IV. Immunoblotting and Detection (all steps are at room temperature, with agitation)

1. Following the electroblotting step, immerse the PVDF membrane in 100% Methanol for 15 seconds, and then allow it to dry at room temperature for 5 minutes.

Note: If Nitrocellulose is used instead of PVDF, this step should be skipped.

2. Block the membrane with 5% non-fat dry milk in TBST for 1 hr at room temperature with constant agitation.

Incubate the membrane with Anti-Cdc42 Antibody, freshly diluted 1:200 to 1:1000 in 5% non-fat dry milk/TBST, for 1-2 hr at room temperature with constant agitation.

Note: To conserve antibody, incubations should be performed in a plastic bag.

3. Wash the blotted membrane three times with TBST, 5 minutes each time.
4. Incubate the membrane with a secondary antibody (e.g. Goat Anti-Mouse IgG, HRP-conjugate), freshly diluted in 5% non-fat dry milk/TBST, for 1 hr at room temperature with constant agitation.
5. Wash the blotted membrane three times with TBST, 5 minutes each time.
6. Use the detection method of your choice. We recommend enhanced chemiluminescence reagents from Pierce.

Example of Results

The following figure demonstrates typical results seen with Cell Biolabs' Cdc42 Activation Assay Kit. One should use the data below for reference only.




Figure 2: Cdc42 Activation Assay. *Left Image:* Cdc42 Immunoblot Positive Control. *Right Image:* Demonstrates Anti-Cdc42 monoclonal antibody specificity by dot blot.




Figure 3: Cdc42 Activation Assay. *Lane 1:* MW Standard. *Lane 2:* 293 cell lysate loaded with GDP and incubated with PAK PBD Agarose beads. *Lane 3:* 293 cell lysate loaded with GTP γ S and incubated with PAK-1 PBD Agarose beads.

References

1. Raftopoulou M., and Hall A. (2004) *Dev Biol.* **265**: 23-32.
2. Bar-Sagi D., and Hall A. (2000) *Cell* **103**: 227-38.
3. Benard, V., Bohl, B. P., and Bokoch, G. M. (1999) *J. Biol. Chem.* **274**, 13198-13204.

Recent Product Citations

1. Li, Q. et al. (2023). Recombinant Slit2 suppresses neuroinflammation and Cdc42-mediated brain infiltration of peripheral immune cells via Robo1-srGAP1 pathway in a rat model of germinal matrix hemorrhage. *J Neuroinflammation.* **20**(1):249. doi: 10.1186/s12974-023-02935-2.

2. Kim, K.B. et al. (2022). WNT5A-RHOA signaling is a driver of tumorigenesis and represents a therapeutically actionable vulnerability in small cell lung cancer. *Cancer Res.* doi: 10.1158/0008-5472.CAN-22-1170.
3. Bijata, M. et al. (2022). Activation of the 5-HT7 receptor and MMP-9 signaling module in the hippocampal CA1 region is necessary for the development of depressive-like behavior. *Cell Rep.* **38**(11):110532. doi: 10.1016/j.celrep.2022.110532.
4. Yuan, Z. et al. (2020). RAB5A promotes the formation of filopodia in pancreatic cancer cells via the activation of cdc42 and β 1-integrin. *Biochem Biophys Res Commun.* **535**:54-59. doi: 10.1016/j.bbrc.2020.12.022.
5. McGuire, S. et al. (2019). Inhibition of fascin in cancer and stromal cells blocks ovarian cancer metastasis. *Gynecol Oncol.* pii: S0090-8258(19)30059-9. doi: 10.1016/j.ygyno.2019.01.020.
6. Hu, H.F. et al. (2018). Comparative Proteomics Analysis Identifies Cdc42-Cdc42BPA Signaling as Prognostic Biomarker and Therapeutic Target for Colon Cancer Invasion. *J Proteome Res.* **17**(1):265-275. doi: 10.1021/acs.jproteome.7b00550.
7. Chang, F. et al. (2017). MTA promotes chemotaxis and chemokinesis of immune cells through distinct calcium-sensing receptor signaling pathways. *Biomaterials.* **150**:14-24. doi: 10.1016/j.biomaterials.2017.10.009.
8. Hong, J. H. et al. (2016). Regulation of the actin cytoskeleton by the Ndel1-Tara complex is critical for cell migration. *Sci Rep.* doi:10.1038/srep31827.
9. Aldinucci, A. et al. (2016). Histamine regulates actin cytoskeleton in human Toll like receptor 4 activated monocyte derived dendritic cells tuning CD4+ T lymphocyte response. *J Biol Chem.* doi:10.1074/jbc.M116.720680.
10. Sherchan, P. et al. (2015). Recombinant Slit2 attenuates neuroinflammation after surgical brain injury by inhibiting peripheral immune cell infiltration via Robo1-srGAP1 pathway in a rat model. *Neurobiol Dis.* **85**:164-173.
11. Sahu, M. et al. (2014). Lens specific RLIP76 transgenic mice show a phenotype similar to microphthalmia. *Exp Eye Res.* **118**:125-134.
12. El-Sayed, F.G. et al. (2014). P2Y₂ nucleotide receptor activation enhances the aggregation and self-organization of dispersed salivary epithelial cells. *Am J Physiol Cell Physiol.* **307**: C83-C96.
13. He, S. et al. (2013). SRGAP1 is a candidate gene for papillary thyroid carcinoma susceptibility. *J. Clin. Endocrinol. Metab.* **98**:E973-E980.
14. Zhang, S. et al. (2008). The tumor suppressor LKB1 regulates lung cancer cell polarity by mediating Cdc42 recruitment and activity. *Cancer Res.* **68**:740-748.
15. Somerville, T. et al. (2006). Identification and characterization of leukemia stem cells in murine MLL-AF9 acute myeloid leukemia. *Cancer Cell* **10**:257-268.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

Contact Information

Cell Biolabs, Inc.
5628 Copley Drive
San Diego, CA 92111
Worldwide: +1 858 271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2006-2026: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.