

Hydroxyproline Assay Kit

Catalog # 6017

For Research Use Only - Not Human or Therapeutic Use

PRODUCT SPECIFICATIONS

DESCRIPTION:	Assay kit to quantify total collagen content (regardless of type or species)
FORMAT:	96-well ELISA Plate with removeable strips
ASSAY TYPE:	Colorimetric assay
ASSAY TIME:	1 hour
STANDARD RANGE:	400 µg/ml to 6.3 µg/ml
NUMBER OF SAMPLES:	Up to 40 (duplicate) colorless samples/plate and up to 20 (duplicate) colored samples/plate
SAMPLE TYPES:	Tissue homogenate and cell culture media
RECOMMENDED SAMPLE DILUTIONS:	Not required
CHROMOGEN:	N/A (read at 530-560 nm)
STORAGE:	-20°C
VALIDATION DATA:	Intra-assay (3.4-6.3%)/Inter-assay (5.1-7.6%)/Spiking Test (91-95%)
NOTES:	Select samples may require a 24-hour hydrolysis step

Hydroxyproline Assay Kit

Catalog # 6017

For Research Use Only - Not Human or Therapeutic Use

INTRODUCTION

Collagen is the major structural protein of the extracellular matrix in many tissues. Hydroxyproline, a major component of collagen, makes up about 13.5% of its amino acid composition. Due to its highly restricted distribution in collagen, the hydroxyproline content accurately reflects the amount of collagen. Therefore, quantitating hydroxyproline has been utilized for evaluating tissue fibrosis or collagen deposition (1 - 3). However, conventional hydroxyproline assays are not useful because they require cumbersome procedures and special tools.


Chondrex, Inc. offers a hydroxyproline assay kit (Cat # 6017) which employs an improved assay system that can be operated with ease and precision using 96-well plates. This kit can quantify the total collagen content in any tissue specimen or tissue homogenate, regardless of collagen type or species. Chondrex, Inc. also provides additional collagen detection kits for different purposes. Firstly, the Sirius Red Total Collagen Detection Kit (Cat # 9062 or 9062P) can be used to quantify solubilized collagen in samples. Secondly, the collagen detection ELISA kits specific to type I or II collagen of various species are used for distinguishing collagen types and species in samples. Lastly, the Sirius Red/Fast Green Collagen Staining Kit (Cat # 9046) is a semi-quantitative assay to determine the total collagen content in tissue sections. For more information, please visit www.chondrex.com or contact support@chondrex.com.

KIT COMPONENTS

Item	Quantity	Amount	Storage
Hydroxyproline Standard (60171)	1 vial	4 mg/ml x 0.5 ml	-20°C
Solution A- Chloramine T Dilution Buffer (60172)	1 bottle	10 ml	-20°C
Solution B - DMAB Dilution Buffer (60173)	1 vial	5 ml	-20°C
10X Chloramine T Concentrate (60174)	1 vial	1 ml	-20°C
2X DMAB (dimethylaminobenzaldehyde) Concentrate (60175)	1 vial	5 ml	-20°C
ELISA Plate	1 each	96-well (8-well strips x 12)	-20°C

Not included: concentrated HCl (10N) and glass screw-thread vials (1-2 ml) with Teflon caps (Example: National Scientific B7999-1)


ASSAY OUTLINE




© 2020 Chondrex, Inc. All Rights Reserved, 6017 7.0

PLATE LAYOUT

Colorless (clear) samples



Colored (turbid) samples



NOTES BEFORE USING ASSAY

NOTE 1: It is recommended that the standard and samples be run in duplicate.

NOTE 2: For partial reagent use, please see the assay protocol's corresponding step for the appropriate dilution ratio. For example, if the protocol dilutes 50 μ l of a stock solution in 10 ml of buffer for 12 strips, then for 6 strips, dilute 25 μ l of the stock solution in 5 ml of buffer. Partially used stock reagents may be kept in their original vials and stored at -20°C for use in a future assay.

NOTE 3: If there are precipitates in the bottles/vials, it is necessary to warm up the bottles/vials in warm water until they completely dissolve.

SAMPLE HYDROLYSIS

1. Weigh 10 mg of a tissue sample in a glass screw-thread vial.
2. Add 100 μ l of distilled water.
3. Mash the tissue sample with a small spatula.

NOTE: 100 μ l of a sample homogenate can be used. Skip steps 1-3 and directly add 100 μ l of the sample homogenate to the vial.

4. Add 100 μ l of concentrated HCl (10N) and tightly screw on the teflon cap.
5. Incubate at 120°C for 24 hours. Mix the sample periodically during incubation.


NOTE: A heat block or dry bath can be used for incubation.

6. Cool down. **Do not open the cap before cooling down.**
7. Use the supernatant for the assay.

NOTE: Black residue is occasionally produced from tissue samples in the hydrolysis process. However, sample hydrolyzation should be complete by the end of the incubation period. If hydrolyzed black residue is still present in the sample, transfer to a microcentrifuge tube and spin at 10,000 rpm for 3 minutes. Collect supernatant.

ASSAY PROCEDURE

- Prepare Standard Dilutions:** Take 100 μl of the 4 mg/ml Hydroxyproline Standard and add to 900 μl of distilled water to make a 400 $\mu\text{g/ml}$ HP standard solution. Then serially dilute it with distilled water. For example, mix 500 μl of the standard (400 $\mu\text{g/ml}$) with an equal volume of distilled water to make a 200 $\mu\text{g/ml}$ solution, and then repeat it five more times to make 100, 50, 25, 12.5, and 6.3 $\mu\text{g/ml}$ standard solutions.



- Prepare Sample Dilutions:** The hydrolyzed samples can be used undiluted. If necessary, the samples can be diluted with 5N HCl. If the samples have color (are not clear), Sample Blank wells should be prepared due to the potential for high background color. See steps 4 and 5 for this process.
- Prepare Chloramine T solution:** Mix 10 μl of the 10X Chloramine T solution with 90 μl of Solution A for each well. For example, 10 samples, 7-point standard, one blank (all in duplicate) will require 3.6 ml of the 1X Chloramine T solution. Mix 360 μl of the 10X Chloramine T solution with 3.24 ml of Solution A.

NOTE: Prepare the 1X solution just before use. Do not store and reuse the mixed solution for the next assay.

- Add Standards and Samples:** Choose 4-1 or 4-2 depending on your samples.

4-1. Colorless Samples: Add 10 μl of standards, distilled water (blank, B) into the purple wells, and samples into the orange wells in duplicate. For example, add 10 μl of sample 1 into the S1 wells, then add 10 μl of sample 2 into the S2 wells. Proceed to Step 5-1.

4-2. Colored Samples: Add 10 μl of standards, distilled water (blank, B) into the purple wells, and samples into the orange and green wells in duplicate. For example, add 10 μl of sample 1 into the S1 and SB1 wells, then add 10 μl of sample 2 into the S2 and SB2 wells. Proceed to Step 5-2.

- Add 1X Chloramine T Solution:**

5-1. Colorless Samples: Add 100 μl of the 1X Chloramine T solution into all wells. Incubate at room temperature for 20 minutes.

5-2. Colored Samples: Add 100 μl of the 1X Chloramine T solution into the purple and orange wells and add 100 μl of Solution A into the green wells. Incubate at room temperature for 20 minutes.

- Prepare DMAB solution:** Mix 50 μl of 2X DMAB solution with 50 μl of Solution B for each well. For example, 10 samples, 7-point standard, one blank (all in duplicate) will require 3.6 ml of the 1X DMAB solution. Mix 1.8 ml of the 2X DMAB solution with 1.8 ml of Solution B.

NOTE: Prepare the 1X solution just before use. Do not store and reuse the mixed solution for the next assay.

- Add 1X DMAB solution:** Add 100 μl of 1X DMAB solution into all wells and incubate at 60°C for 30 minutes.

NOTE: Mix the plate by tapping gently or using a plate shaker.

- Read Plate:** Read the OD values at 530-560 nm. If the OD values of samples are greater than the OD values of the highest standard, re-assay the hydrolyzed samples at a higher dilution using 5N HCl.

CALCULATING RESULTS

1. Average the duplicate OD values for the blank, standards, test samples, and sample blanks (if used).
2. Subtract the "blank" (B) values from the averaged OD values in step 1. If sample blanks were used, also subtract the averaged sample blank values from their corresponding averaged sample values.
3. Plot the OD values of standards against the concentration of the standards ($\mu\text{g/ml}$). Using a log/log plot will linearize the data. Figure 1 shows a representative experiment where the standard range is 6.3 - 400 $\mu\text{g/ml}$.
4. The concentration of Hydroxyproline ($\mu\text{g/ml}$) in the hydrolyzed samples can be calculated using regression analysis. Multiply the results by the dilution factors if the hydrolyzed samples were diluted.
5. Choose one of the following equations depending on the sample type:

5-1 Solid Samples: Hydroxyproline level in a tissue sample ($\mu\text{g/mg}$) is determined by the following equation:

$$\frac{\text{Hydroxyproline } (\mu\text{g/ml}) \times (\text{Distilled Water Volume (ml)} + \text{HCl Volume (ml)})}{\text{Sample Weight (mg)}} = \text{Hydroxyproline level } (\mu\text{g/mg}) \text{ in the sample}$$

5-2 Solution or Homogenate: Hydroxyproline level in a solution sample (mg/ml) is determined by the following equation:


$$\frac{\text{Hydroxyproline } (\mu\text{g/ml}) \times (\text{Sample Volume (ml)} + \text{HCl Volume (ml)})}{\text{Sample Volume (ml)}} = \text{Hydroxyproline level } (\mu\text{g/ml}) \text{ in the sample}$$

6. Hydroxyproline levels can be converted into collagen levels with the following equation (4):

$$\text{Hydroxyproline level } (\mu\text{g/mg or } \mu\text{g/ml}) \times \frac{100}{13.5} = \text{Collagen level } (\mu\text{g/mg or } \mu\text{g/ml})$$

Note: Hydroxyproline accounts for 13.5% of the collagen amino acid composition.

Figure 1 - A Typical Standard Curve for the Hydroxyproline Assay Kit



ASSAY VALIDATION

Table 1 - Reproducibility Data for the Hydroxyproline Assay Kit

Test	Mouse Kidney	Hydroxyproline 200 µg/ml	Hydroxyproline 12.5 µg/ml
Intra-Assay CV (%)	5.1	6.3	3.4
Inter-Assay CV (%)	5.1	6.1	7.6
Spike Test* (%)	91-95%	-	-

*Standard was mixed with known amounts of mouse kidney samples or hydroxyproline solution

TROUBLESHOOTING

For frequently asked questions about assays and ELISAs, please see Chondrex, Inc.'s [Assay FAQ](#) for more information.

REFERENCES

1. N. Blumenkrantz, G. Asboe-Hansen, A Quick and Specific Assay for Hydroxyproline. *Anal Biochem* **55**, 288-91 (1973).
2. G. Reddy, C. Enwemeka, A Simplified Method for the Analysis of Hydroxyproline in Biological Tissues. *Clin Biochem* **29**, 225-9 (1996).
3. C. Rogers, J. Kimmel, M. Hutchin, H. Harper, A Hydroxyproline Method of Analysis for a Modified Gelatin in Plasma and Urine. *J Biol Chem* **206**, 553-9 (1954).
4. R. Neuman, M. Logan, The Determination of Hydroxyproline. *J Biol Chem* **184**, 299-306 (1950).